



The 7000C Dynamic Headspace module is a plug and play option for CDS 7000C concentrator. This module, when installed with 7000C concentrator, adds dynamic headspace capability to PAL RTC/RSI rail. The setup is capable of performing dynamic headspace analysis in a closed-system, which is USP 467, EPA 5035A and 8260C soil compliant.



User Interface of PAL DHS Plugin



- High Temp Isothermal Zones**
1. Max 300 °C Sample Vial Heater
 2. Max 300 °C Valve Oven
 3. Max 350 °C Transferline

- Seamless Upgrade**
1. Plug and Play Module
 2. Easy Mount on PAL System
 3. TCP/IP Based

- Build-in Magnetic Stirrer**
1. Select from Three Speeds
 2. 8260C Compliant in Soil Mode

